



Be Right™



Test strips pH , 4-9 pH units

Kat. číslo produktu:	2745650
CZK Cena bez DPH:	Kontaktujte nás
Dostupné	

Most natural waters range from pH 4 to pH 9, but commonly are above pH 7 because of carbonates and bicarbonates (alkalinity). The pH of water used in industry, boiler feed water, and swimming pools usually is kept within a narrow range. Contains 50 tests.

Easy to use, disposable, and inexpensive

Great for quick evaluations of pH levels in the lab or in the field

Technické údaje

Jednotky:	50 tests
Nejmenší dílky přírůstku:	Steps: 4, 5, 6, 7, 8, 9 pH
Parametr:	pH, mid range
Platforma:	Test Strip
Rozsah měření:	4 - 9 pH units
Typ testovací soupravy:	Test Strips