



Soxhlet thimbles, cellulose, 33 x 80 mm, 25 pk

Kat. číslo produktu:

2966700

ZASTARALÁ POLOŽKA

Tato položka již není dostupná.

Xenosep[®] extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep[®] extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

Technické údaje

Materiál:	Cellulose
Obsah balení:	25 /pk
Velikost póru:	10 µm