



BAG, WHIRL-PAK* 5-1/2X9 PK/10

Kat. číslo produktu:	1437297
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection.

Technické údaje

Kapacita:	720 mL
Obsah balení:	10