




Be Right™



SAC950 Sample changer for automation of large sample series (Radiometer Analytical)

Kat. číslo produktu: R41T044

ZASTARALÁ POLOŽKA
Tato položka již není dostupná.

Effortless sample handling due to total automation

Extended sample changer for automation of large sample series, multi-radius sample trays for up to 70 samples

Compact and adaptable

Flexible and programmable

Quick and efficient

Safe and clean

Technické údaje

Délka:	686 mm
Hmotnost:	27 kg
Rozměry (V x Š x H):	70 mm x 580 mm x 686 mm
Typ:	Autosampler