



Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Kat. číslo produktu:	6788302
CZK Cena bez DPH:	Kontaktujte nás
Dostupné	

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Technické údaje

Jaký je obsah balení:	Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions
Product Kit:	5500 sc Silica Analyser, 2/4 Channels

Jaký je obsah balení

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions