



## PHC3085-8 kombinovaná pH elektroda, teplotní čidlo

Kat. číslo produktu:

E16M342

### **ZASTARALÁ POLOŽKA**

Tato položka již není dostupná.

### **The right probe for your application.**

To meet your analysis needs, Radiometer Analytical electrodes are available with a variety of lengths, diameters and reference systems. There is a Radiometer electrode for every application.

### **Glass sensors made for long-life and durable use.**

Radiometer Analytical glass pH electrodes are made with highly shock-resistant glass, making them extremely tough. The average lifetime of a correctly maintained Radiometer glass pH electrode is approximately two years.

---

## **Technické údaje**

|                |             |
|----------------|-------------|
| Materiál:      | Body: epoxy |
| Typ elektrody: | pHC3085-8   |
| Záruka:        | 6 měsíců    |