



Iron color disc test kit, model IR-18B

Kat. číslo produktu:

146401

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dostupné



Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Color disc/ 1,10 Phenanthroline
Parameter:	Iron
Parametr:	Iron, medium range - as Fe**
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	IR-18B
Rozsah měření:	High Single Range
Smallest increment:	0.25
Standardní rozsah měření:	0.25 - 7 mg/L Fe
Styl kufříku:	B
Typ:	S jednou stupnicí
Typ testovací soupravy:	IR-18B