



Total Hardness Test Kit, Model 5-EP

Kat. číslo produktu:

145400

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů



Water hardness is caused almost entirely by calcium and magnesium ion. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With the Model 5-EP test kit, simply count the number of drops required to change the solution color to determine hardness concentration. Utilizes drop count titration with a range of 1-30 gpg. Contains approximately 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Drop count titration/ EDTA-powder pillows
Parametr:	Hardness, Total - As CaCO_3
Platforma:	Drop Count Titrator
Počet testů:	100
Rozsah měření:	1 - 30 gpg CaCO_3
Smallest increment:	1 gpg
Styl kufříku:	B
Typ testovací soupravy:	5-EP