



Hardness (Total & Calcium) Test Kit, Model HAC-DT

Kat. číslo produktu:

2063900

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Technické údaje

Jaký je obsah balení:	Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.
Metoda:	Digital Titrator/EDTA
Parametr:	Hardness, total & Calcium - As CaCO_3
Platforma:	Digital Titrator
Počet testů:	100 Ca +100 T
Rozsah měření:	10 - 4,000 mg/L CaCO_3
Smallest increment:	0.1-10
Styl kufříku:	F
Typ testovací soupravy:	HAC-DT

Jaký je obsah balení

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.