



Test kit, free chlorine, model CN-70F

Kat. číslo produktu:

1454201

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Quickly and accurately measure free chlorine levels in the field or in the lab

Technické údaje

Metoda: Color disc/DPD

Parameter: Chlorine Free

Parametr: Chlorine, low range - As free Cl₂

Platforma: Barevný disk

Počet testů: 200

Rozsah měření: Dual Range

Rozsah měření (2): 0 - 3,4 mg/L Cl₂

Smallest increment: 0.02

0.1

Styl kufříku: D

Typ testovací soupravy: CN-70F²

Vysvětlivka: *mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.