



Dissolved oxygen test kit, model OX-2P

Kat. číslo produktu:

146900

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Nebezpečné



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

Utilizes drop count titration - a simple, economical titration method

Technické údaje

Metoda:	Drop count titration/ Modified Winkler
Parametr:	Dissolved Oxygen - As DO
Platforma:	Drop Count Titration
Počet testů:	100
Rozsah měření:	0.2 - 4.0 mg/L & 1 - 20 mg/L O ₂
Smallest increment:	0.2
Styl kufříku:	D
Typ testovací soupravy:	OX-2P