



Be Right™



Bottle, Wash, Nalgene®, Narrow Mouth, 500 mL, Red Cap/ Stem, 6/pk

Kat. číslo produktu:

2928204

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Technické údaje

Barva:	Čirá
Jednotka:	6 pk
Kapacita:	500 mL
Materiál:	LDPE
Obsah balení:	6 pk
Stem:	Red Cap/Stem