



BEAKER, PP LOW FORM 100ML

Kat. číslo produktu:

108042

CZK Cena bez DPH:

Kontaktujte nás

Dostupné

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Technické údaje

Jednotky:	1 each
Kapacita:	100 mL
Materiál:	Polypropylene, Low Form
Materiál skříně:	Polypropylene
Obsah balení:	1 each
Průměr:	63.0 mm with spout