



pHC2401-7 Combined pH electrode, Red Rod, ring diaphragm, Type 7

Kat. číslo produktu:

E16M336

ZASTARALÁ POLOŽKA

Tato položka již není dostupná.

Navrhované náhrady:

- E16M400
- 945-252

Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: annular ring. With plug Type 7.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction