



DRB 200-2 Dry thermostat with 2 blocks, 15x16 mm, 15x16 mm

Kat. číslo produktu: LTV082.99.44001

CZK Cena bez DPH: Kontaktujte nás
Dostupné

Thermostat DRB 200 – the simple solution for all digestions

Pre-programmed for standard digestions (e.g. for P_{tot} , TN or metals) and freely programmable for user specific digestions.

Equipped with a large readable display and a digital timer with an automatic switch-off and acoustic signal. Operational safety guaranteed by protection lid covering the heating block, by the insulated external shell and by an integrated anti-overheating feature. Cuvettes and reaction vessels can be digested at large scale of different temperatures and time settings. Temperature stability and reproducibility ensured.

Great flexibility

Excellent reproducibility

Simple to use

Technické údaje

Bezpečnostní kontroly:	CE, GS, cTUVus
Displej:	5 x 7 dot matrix
Heating system:	2 Heating Blocks
Hmotnost:	2.8 kg
Jazyky manuálu:	English, French, German, Italian, Spanish, Dutch, Swedish, Danish, Polish
Max. provozní vlhkost:	90 %
Počet kyvet:	15 x 16 mm diameter
Počet kyvet (2):	15 x 16 mm diameter
Podmínky skladování:	-40 - 60 °C
Požadavky na napájení (V):	115/220 V AC
Programy ohřevu:	Preprogrammed for 40°C, 100°C, 150°C, 165°C and freely selectable 37-165°C, 1-480 min.
Provozní teplotní rozsah:	10 - 45 °C
Rychlost ohřevu:	From 20 - 150 °C in 10 minutes
Síťové napájení:	115 V - 230 V +5%/-15%
Teplotní stabilita:	± 1 °C v souladu s metodami EN, ISO, EPA
Typ displeje:	LCD
Uživatelské programy:	6 free temperature/time
Uživatelské rozhraní:	English, French, German, Italian, Spanish, Dutch, Swedish, Polish

Velikost displeje:	2 x 16 character
Záruka:	24 měsíců