



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Kat. číslo produktu:

108076

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Technické údaje

Jednotky:	6 /pk
Kapacita:	250 mL
Materiál:	Polypropylene, Low Form
Materiál skříně:	Polypropylene
Průměr:	83.6 mm with spout