



Be Right™



Flask, Volumetric, Polypropylene, 50 mL

Kat. číslo produktu:

1406041

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Technické údaje

Certifikace:	Not Applicable
Jednotky:	1 each
Kapacita:	50 mL
Materiál skříně:	Polypropylene