



Test kit, copper (free) test color disc, model CU-5

Kat. číslo produktu:

1421300

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Copper may exist in natural waters and effluents as a soluble salt or as suspended solids. A small amount is essential for plants and animals, and concentrations exceeding 0.1 mg/L are useful for controlling algae and plankton growth. Quantities ranging from 0.02 to 0.1 mg/L are toxic for some fish, so its use for treating fish disease requires careful monitoring. This kit is USEPA-accepted for waste water and potable water and approved by some states for potable water under specific conditions (Check with your regulatory agency). Range: 0.1-5 mg/L; contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Color-disc/Bicinchoninate
Parameter:	Copper (free)
Parametr:	Copper - As free Cu
Platforma:	Barevný disk
Počet testů:	100 F
Přístroj:	CU-5
Rozsah měření:	0 - 4 mg/L Cu
Smallest increment:	0.1
Standardní rozsah měření:	0 - 4 mg/L Cu
Styl kufříku:	B
Typ testovací soupravy:	CU-5