



Be Right™



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 250 mL, pk/12

Kat. číslo produktu:	1472456
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Polyethylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Technické údaje

Body Material:	Plastic
Kapacita:	250 mL
Materiál:	PP
Obsah balení:	12 /pk
Package size UOM:	12pk