



**Be Right™**



## Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 500 mL, 12/pk

**Kat. číslo produktu:** 1472459

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

---

### Technické údaje

|                   |         |
|-------------------|---------|
| Body Material:    | Plastic |
| Kapacita:         | 500 mL  |
| Materiál:         | PP      |
| Obsah balení:     | 12 /pk  |
| Package size UOM: | 12pk    |