



Hydrazine color disc test kit, model HY-2

Kat. číslo produktu:

184900

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method with a range of 0.02-1 mg/L, this kit contains 300 tests.

Measures Hydrazine, a common oxygen scavenger in boiler water, using the simple color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Color disc/ p-Dimethylamino- benzaldehyde
Parameter:	Detergents
Parametr:	Hydrazine - As N_2H_4
Platforma:	Barevný disk
Počet testů:	300
Přístroj:	HY-2
Rozsah měření:	0 - 1 mg/L N_2H_4
Smallest increment:	0
Standardní rozsah měření:	0 - 1 mg/L N_2H_4
Styl kufříku:	D
Typ testovací soupravy:	HY-2