



Test kit, hardness (total), model HA-DT

Kat. číslo produktu:

2063600

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Utilizing the Hach Digital Titrator, the Model HA-DT kit is used to determine hardness concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Technické údaje

Metoda:	Digital Titrator/EDTA
Nejmenší dílky přírůstku:	0.1-10
Parametr:	Hardness, total - As CaCO_3
Platforma:	Digital Titrator
Počet testů:	100
Rozsah měření:	10 - 4,000 CaCO_3
Styl kufříku:	F
Typ testovací soupravy:	HA-DT