



**Be Right™**



## Test kit, chloride, model CD-51

**Kat. číslo produktu:**

**208601**

CZK Cena bez DPH:  
Datum odeslání neuvedeno

Kontaktujte nás

**Nebezpečné**



Found in nearly all natural waters, chloride affects human taste above 250 mg/L. High levels inhibit plant growth, and many industrial processes need to limit chloride concentration. The model CD-DT, multi-range kit, is supplied with the demineralizer bottle for high-range sample dilution. Contains 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Technické údaje

Metoda:	Drop count titration/
	Silver nitrate
Parametr:	Chloride - As Cl <sup>-</sup>
Platforma:	Drop Count Titration
Počet testů:	100
Rozsah měření:	500 - 10000 mg/L & 5000 - 100000 mg/L Cl <sup>-</sup>
Smallest increment:	500
	5000
Styl kufříku:	B
Typ testovací soupravy:	CD-51
Vysvětlivka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.