



Be Right™



Bottle, Sampling, Low Density, Polyethylene, 250 mL, 12/pk

Kat. číslo produktu:

2087076

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Technické údaje

Body Material:	Low Density Plastic
Kapacita:	250 mL
Materiál:	LDPE
Obsah balení:	12 /pk
Package size UOM:	12pk