



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 500 mL

Kat. číslo produktu:	2089679
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Technické údaje

Jednotky:	6/pk
Kapacita:	500 mL
Materiál:	Glass Screw Cap