



## Flask, Erlenmeyer, glass w/screw cap, 6/pk, 1000 mL

Kat. číslo produktu: 2089683

**ZASTARALÁ POLOŽKA**  
Tato položka již není dostupná.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

---

### Technické údaje

|           |                 |
|-----------|-----------------|
| Jednotky: | 6/pk            |
| Kapacita: | 1000 mL         |
| Materiál: | Glass Screw Cap |