



Be Right™



Flask, Erlenmeyer, Polycarbonate Capacity 125 mL

Kat. číslo produktu:	2089743
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Graduated (Molded-in graduations) flask features: clarity, toughness and strength. Autoclavable; however, repeated autoclaving may result in loss of mechanical strength. #5 Neoprene stoppers (P/N 14808-05; pk/12) available separately. 125 mL.

Technické údaje

Jednotky:	1 each
Kapacita:	125 mL
Materiál:	Polycarbonate