



Test kit, formaldehyde, model FM-1

Kat. číslo produktu:

2183100

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Nebezpečné



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Drop count titration/ Thymolphthalein
Parametr:	Formaldehyde - As CH_2O
Platforma:	Drop Count Titration
Počet testů:	100
Rozsah měření:	0,05 - 1% & 0,5 - 10%
Smallest increment:	0.05% 0.5%
Styl kufříku:	D
Typ testovací soupravy:	FM-1