



Reagent set, LR molybdenum, 5 mL

Kat. číslo produktu:

2439500

Nebezpečné

ZASTARALÁ POLOŽKA

Tato položka již není dostupná.



For determination of Low Range Molybdate Molybdenum by the ternary complex method, using visual test kits. Replacement reagent set for MO-LR Test Kit (2359300). Range: 0.25 - 3 mg/L Mo. Sample Size: 5 mL. Set contains: Molybdenum 1 Reagent Powder Pillows and Molybdenum 2 Reagent Solution for approximately 100 tests.

Technické údaje

Metoda:	Ternary Complex
Rozsah měření:	0.25 - 3 mg/L Mo
Stojan na vzorek:	5 mL