



Be Right™



Flask, Erlenmeyer, wide mouth, 125 mL, Subdivisions 50 mL

Kat. číslo produktu:	2489443
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included.

Technické údaje

Jednotky:	1 each
Kapacita:	125 mL
Materiál:	Glass Wide Mouth