



## Bottle, sampling, presterilized, round, 12/pk

Kat. číslo produktu: 2495012

**ZASTARALÁ POLOŽKA**  
Tato položka již není dostupná.

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Works well with P-A Broth Ampules and PathoScreen™ Medium.

---

### Technické údaje

|                        |       |
|------------------------|-------|
| Dechlorinating Agent?: | No    |
| Obsah balení:          | 12/pk |