



Be Right™



Sample container, 240 mL

Kat. číslo produktu:	2552411
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Plastic, leak-proof. Corning Snap-Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with unique sealing mechanism so sample won't leak.

Technické údaje

Kapacita:	240 mL
Obsah balení:	140 /pk
Průměr:	80 mm (approx)