



Test kit, chloride, model CDS-DT

Kat. číslo produktu:

2580600

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

Nebezpečné



The single-parameter Model CDS-DT uses the Hach Digital Titrator with Silver Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-10,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Technické údaje

Metoda:	Digital Titrator/ Silver nitrate
Parametr:	Chloride - As Cl^- , 10-10,000 mg/L Cl^-
Platforma:	Digital Titrator
Počet testů:	100
Rozsah měření:	10 - 10,000 mg/L Cl^-
Rozsah měření (2):	to
Smallest increment:	10
Styl kufříku:	F
Typ testovací soupravy:	CDS-DT

Požadované příslušenství

- Pipeta, variabilní objem, 1,0-5,0 mL (Item BBP065)