



Microfil funnel & filter, 100 mL, 0.45 µm 47 mm white gridded (millipore)

Kat. číslo produktu:	2586300
CZK Cena bez DPH:	Kontaktujte nás
Dostupné	

Designed to eliminate cleaning and sanitization of filter funnels between uses. Graduated (20, 50, and 100 mL) 100 mL funnels are sterile and ready-to-use. Funnel creates leak-proof seal when seated onto Microfil Stand Alone Support. Excellent for testing water samples from potable, non-potable, source, ambient, industrial, food, beverage, pharmaceutical and wastewater sources. Microfil funnels come with 150 Millipore filters (white; gridded; cellulose and mixed ester composition; 0.45 µm pore size; 47 mm diameter; sterile; individually packed). Microfil V includes integrated filter for added convenience and can be used with standard No. 8 stopper on any filter holder or manifold.

Technické údaje

Kapacita:	100 mL
Materiál:	Mixed-Cellulose Ester
Obsah balení:	150
Průměr:	47 mm
Velikost póru:	0.45 µm