



## Buret, Class A, Glass, 25 mL

|                          |                 |
|--------------------------|-----------------|
| Kat. číslo produktu:     | 2636540         |
| CZK Cena bez DPH:        | Kontaktujte nás |
| Datum odeslání neuvedeno |                 |

Glass, funnel-fill buret. Certified Class A (ASTM Specification E287) with permanent serial number and certificate for NIST traceability. PTFE stopcock does not require grease. Dust cap included. 10 mL buret meets requirements for potentiometric titrations (ASTM D664).

---

### Technické údaje

|                   |               |
|-------------------|---------------|
| Kapacita:         | 25 mL         |
| Materiál skříně:  | Class A Glass |
| Speciální funkce: | Buret         |