



Test strips total Hardness, 0-425 mg/L, 250 tests, individually wrapped

Kat. číslo produktu: 2793844

ZASTARALÁ POLOŽKA
Tato položka již není dostupná.

Water hardness is caused almost entirely by calcium and magnesium ions. Other di- and trivalent metals have a similar effect, but usually are not present in high enough concentration in potable waters to cause problems. Hardness increases soap consumption in laundries and causes scale in boilers. With Hach Water Quality Test Strips, technicians in the field can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Hach test strips are also used in laboratories all over the world for pre-test screening-to detect the presence of materials that might interfere with lab testing. Strips are individually wrapped for maximum portability and stability. Contains 250 tests.

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Technické údaje

Jednotky:	250 tests (individually wrapped)
Nejmenší dílky přírustku:	Steps: 0, 25, 50, 120, 250, 425 ppm
Parametr:	Hardness, Total (as CaCO ₃)
Platforma:	Test Strip
Počet testů:	250
Rozsah měření:	0 - 425 ppm
Rozsah měření (2):	0 - 25 gpg
Smallest increment:	Intervalos de: 0, 1.5, 3, 7, 15, 25 gpg
Typ testovací soupravy:	Test Strips
Vysvětlivka:	*mg/L unless otherwise noted; ppb = µg/L; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.