



Test kit color disc chlorine dioxide, model CLD-2

Kat. číslo produktu:

2842800

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů



Monitor ClO_2 levels on-site with the NEW Chlorine Dioxide Test Kit. The Chlorine Dioxide Test Kit uses the reliable and EPA-accepted DPD method to monitor ClO_2 concentration, while eliminating chlorine interference through the addition of glycine. With 100 tests per kit and available range of 0-2 mg/L and resolution of 0.1 mg/L, operators can accurately monitor around the 0.8 mg/L regulated limit.

Simple, cost-effective method for maintaining optimum chlorine dioxide levels in your plant

USEPA-accepted for waste water and potable water

Approved by some states for potable water under specific conditions (check with your regulatory agency)

Technické údaje

Metoda:	Color Disc/DPD
Parameter:	Chlorine Dioxide
Parametr:	Chlorine Dioxide
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	CLD-2
Rozsah měření:	0 - 2,0 ClO_2
Smallest increment:	0.1
Standardní rozsah měření:	0 - 2.0 ppm ClO_2
Styl kufříku:	D
Typ testovací soupravy:	CLD-2 ²