



Membrane filters millipore S-pak type HA, 47 mm, 0.45 um pore size, pk/600

Kat. číslo produktu: 2936500
CZK Cena bez DPH: Kontaktujte nás
Datum odeslání neuvedeno

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 um. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Technické údaje

Materiál:	Mixed-Cellulose Ester
Obsah balení:	600
Průměr:	47 mm mm
Velikost póru:	0.45 µm