



QL filters, pore size 8 - 12 μm , diameter 90 mm, circular, 100/pk, qualitative

Kat. číslo produktu: 50655

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Technické údaje

Materiál:	Cellulose
Obsah balení:	100 /pk
Průměr:	90 mm