



QL filters, pore size 8 - 12 μm , diameter 125 mm, circular, 100/pk, qualitative

Kat. číslo produktu:	50657
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Technické údaje

Materiál:	Cellulose
Obsah balení:	100 /pk
Průměr:	125 mm