



Be Right™



Funnel, separatory, glass, 500 mL

Kat. číslo produktu:	52049
CZK Cena bez DPH:	Kontaktujte nás
Datum odeslání neuvedeno	

Pear-shaped (Squibb) separatory funnel with PTFE stopcock. Clear glass body allows for easy detection of layers between immiscible liquids, all the way to the stopcock. Highly chemical resistant and fully autoclavable. High density linear polyethylene s

Technické údaje

Kapacita:	500 mL
Obsah balení:	1 each