



Glass fiber filter, pore size 1.5 µm, 47 mm diameter, 100 pieces

Kat. číslo produktu:	253000
CZK Cena bez DPH:	Kontaktujte nás
Dostupné	

These filters are recommended for use in Standard Methods 2540 C, 2540 D, and 2540 E.

Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results comparable to 934-AH® brand filters. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater. These filters can withstand temperatures over 500°C which is an especially important feature for laboratories performing Fixed and Volatile Solids analysis of water and wastewater. ← 934-AH is a registered trademark of Whatman, Inc./GE Healthcare.

Technické údaje

Materiál:	Glass Microfiber
Obsah balení:	100 /pk
Průměr:	47 mm
Velikost póru:	1.5 µm