



LABCOM Easy Komunikační software, pro sensION+ SLP, CD, kabel, adaptér USB

Kat. číslo produktu: LZW8997.99

ZASTARALÁ POLOŽKA
Tato položka již není dostupná.

LabCom Easy software, to obtain data of Sension+ GLP instruments from a PC. Kit including CD and cable with 9-contact connector and USB adapter.

- On line graph view of measured value
- Results table personalisation, the user selects the displayed information
- Storage of calibration and measuring results in a data base where the data can be viewed or printed at any moment
- Print report of individual sample or results summary for several samples

Data export to Excel

Convenient Data Management.

Labcom easy allows the acquisition and the management of the obtained results from Sension+ Laboratory instruments.