



db CDC304 sonda konduktivity

Kat. číslo produktu: 945-002

ZASTARALÁ POLOŽKA
Tato položka již není dostupná.

Vodná a nevodná média.

The right conductivity cell for your application.

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

Technické údaje

Jaký je obsah balení: Includes: Radiometer Analytical CDC304 conductivity cell with Certificate of Conformity.

Jaký je obsah balení

Includes: Radiometer Analytical CDC304 conductivity cell with Certificate of Conformity.