



## Radiometer Analytical ISE25Ca-9 Calcium Ion Selective Electrode (ISE) (half-cell, screw cap)

Kat. číslo produktu: E41M002

**ZASTARALÁ POLOŽKA**  
Tato položka již není dostupná.

The Radiometer ISE25Ca-9 Calcium Ion Selective Electrode (ISE) has a plastic body, a refillable internal electrolyte (article no. S41M001), and calcium ion ( $\text{Ca}^{2+}$ ) specific, replaceable PVC membrane (article no. E91M001). The ISE25Ca-9 is a half-cell ISE and requires an external reference. The Radiometer REF201 (article no. E21M009), REF251 (article no. E21M001) or REF451 (article no. E21M005) reference electrode is recommended. This ISE is a screw cap electrode and requires the CL114 cable (part no. A94L114) for compatibility with Hach Sension+ benchtop meters (or any ion meter with a BNC connection).

---

### Technické údaje

|                           |                |
|---------------------------|----------------|
| Délka:                    | 103 mm         |
| Provozní teplotní rozsah: | 0 - 50 °C      |
| Rozměry (H x D):          | 12 mm x 103 mm |