



Iron color disc test kit, model IR-18A

Kat. číslo produktu:

146500

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Color disc/ 1,10 Phenanthroline
Parameter:	Iron
Parametr:	Iron, medium range - As Fe**
Platforma:	Barevný disk
Počet testů:	100
Přístroj:	IR-18A
Rozsah měření:	Low Single Range
Smallest increment:	0.02
Standardní rozsah měření:	0.02 - 1.00 mg/L Fe
Styl kufříku:	B
Typ:	S jednou stupnicí
Typ testovací soupravy:	IR-18A