



Test Kit, chlorine (free & total) & pH, model CN-67

Kat. číslo produktu:

1411100

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dodání do 2 týdnů



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

One of Hach's most popular test kits

Includes color discs and a comparator for monitoring chlorine and pH

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Aplikace:	Pool and Spa
Parametr:	Chlorine (Free & Total) and pH
Platforma:	Barevný disk
Počet testů:	300
Přístroj:	CN-67
Rozsah měření:	6.6 - 8.4 pH units
Rozsah měření (2):	0 - 3,4 mg/L