



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Kat. číslo produktu:

143601

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Dostupné



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Technické údaje

Metoda:	Drop count titration/ Sodium hydroxide Phenolphthalein
Parametr:	Carbon Dioxide - As CO ₂
Platforma:	Drop Count Titration
Počet testů:	200
Rozsah měření:	2 - 40 mg/L
Rozsah měření (2):	1.25 - 25 mg/L
Rozsah měření (3):	5 - 100 mg/L CO ₂
Smallest increment:	1.25 2 5 D
Styl kufříku:	CA-23
Typ testovací soupravy:	
Vysvětlivka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.