



## Test kit cyanide, model CYN-3

**Kat. číslo produktu:**

**201002**

**CZK Cena bez DPH:**

**Kontaktujte nás**

**Datum odeslání neuvedeno**

**Nebezpečné**



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures.

**Simple, low-cost method using a stepped color comparator**

**Achieve fast results with an easy-to-carry test kit**

---

### Technické údaje

|                           |                                   |
|---------------------------|-----------------------------------|
| Metoda:                   | Color disc/                       |
|                           | Pyridine-Pyrazolone               |
| Parameter:                | Cyanide                           |
| Parametr:                 | Cyanide - As free CN <sup>-</sup> |
| Platforma:                | Barevný disk                      |
| Počet testů:              | 100                               |
| Přístroj:                 | CYN-3                             |
| Rozsah měření:            | 0 - 0,3 mg/L CN <sup>-</sup>      |
| Smallest increment:       | 0.02                              |
| Standardní rozsah měření: | 0 - 0.3 mg/L CN <sup>-</sup>      |
| Styl kufříku:             | R                                 |
| Typ testovací soupravy:   | CYN-3                             |