



Test kit, alkalinity, model AL-DT

Kat. číslo produktu:

2063700

Nebezpečné

CZK Cena bez DPH:

Kontaktujte nás

Datum odeslání neuvedeno



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Determine phenolphthalein and total alkalinity by visual endpoint titration. Using the Digital Titrator, this contains 100 tests with a range of 10-4,000 mg/L.

Digital Titrator and standardized titration cartridges provide laboratory accuracy for field titrations

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Technické údaje

Metoda:	Digital Titrator/Sulfuric acid
Nejmenší dílky přírustku:	0.1-10
Parametr:	Alkalinity - (P) & Total (MO) as CaCO ₃
Platforma:	Digital Titrator
Počet testů:	100
Rozsah měření:	10 - 4,000 mg/L CaCO ₃
Styl kufříku:	F
Typ testovací soupravy:	AL-DT
Vysvětlivka:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.